

Simple Interventions. Extraordinary Outcomes.



COMFORT PERSONAL CLEANSING® WARMER



IMPROVED FEATURES FOR IMPROVED RELIABILITY

Preparing Comfort Personal Cleansing Products for use has never been easier. Our warmers now include improved features such as a digital display for easy viewing of information, utilization tracking, adjustable temperature range and the ability to record audio cues.

- **Keep better track of product inventory**
The warmer's new utilization tracking feature allows you to see how many packages are being used. The warmer automatically records when a package is removed for use.
- **Keep products at just the right temperature**
Adjustable temperature range of 100°F to 125°F helps you warm product to protocol specifications and patient comfort.
- **Alert staff to proper product use and more with audio cues**
Recordable audio cues mean you can record clinical reminders, restocking information, Do Not Flush reminders and more!

**TAKE
FIRST**

NO MORE GUESSING WHICH PRODUCTS ARE WARM AND READY FOR USE.

READY

Each product is warmed to a consistent, reliable temperature. Meanwhile, lighted indicators automatically show which products are warm and ready, and which ones are not. They even promote efficient "first in/first out" stock rotation.

NOT READY

ADDITIONAL FEATURES:

- **Improved product sensing**
Accurately records which products are ready for use and which need to be disposed.
- **Clear view of products and indicators**
Front-loading, see-through doors make it easy to see if warmer is fully stocked and if products are ready for use.
- **Compatible with all warmable products from Sage.**
- **Internationally compatible**
UL/IEC/CE approved – it meets most international standards for safety for medical electrical equipment.
- **Quiet power supply contains no moving parts.**
- **Conduction heating mechanism for more reliable warming of products.**
- **Sleek, modern design.**

CAUTION: It is very important to check the temperature of the washcloths. If there is any indication that washcloths are excessively hot, do not use. Remember, gloves diminish your sensitivity to heat. During bathing, communicate with your patient to make sure each washcloth is at a comfortable temperature. Closely monitor heat level with infants, unresponsive or skin sensitive individuals.



Digital display for easy viewing of information

Alert staff to proper product use and more with audio cues



Comfort Personal Cleansing® Cart

- Convenient, point-of-use work station and storage solution.
- Accommodates one warmer and four dispenser boxes of Comfort Personal Cleansing products.

PRODUCT ORDERING:



28-COUNT WARMER
19.625" (H) x 27.5625" (W) x 15.375 (D)
1/case
Reorder #7938



CART
1/case
Reorder #7920